

Instructions for Use / Datasheet
SwiftFISH Rapid Hybridization Buffer Protocol

Rev Date: Dec 19 2016

Rev: 2

1 of 1

SwiftFISH Rapid Hybridization Buffer Protocol

Notes

- The SwiftFISH Rapid Hybridization Buffer decreases hybridization time from 16 hours (traditional hybridization buffer) to as little as 1 hour.
- Probe intensity is equivalent to use of traditional hybridization buffer
- Works for the following samples types: peripheral blood (PB), bone marrow (BM) and formalin-fixed paraffin-embedded (FFPE)
- Works with all FISH probes: controls, gene specific, custom FISH probes
- Further optimization of the protocol may be required.

Required Reagents & Equipment (Not Supplied)

Hotplate 70%, 85%,

70%, 85%, 100% Ethanol

Humidified Slide Chamber

Wash Solution 1 (WS1) -0.3% Igepal (Sigma CA-630) or NP-40 / 0.4 x SSC

Wash Solution 2 (WS2) – 0.1% Igepal (Sigma CA-630) or NP-40 / 2 x SSC

DAPI with Antifade

Hybridization Setup

- 1. Take buffer from fridge and pre-warm to 37°C.
- 2. Mix buffer by tapping/vortexing, making sure there is no precipitate.
- 3. Place slides on a warm plate (45°C) for 20 minutes.
- 4. Dehydrate the slides in Ethanol series 70%, 85%, and 100%, 2 minutes each. Let dry.
- 5. Add 10 µl probe mixture to slide (2 µl probe + 8 µl SwiftFISH Rapid Hybridization Buffer) using a pipette tip with the end cut.
- 6. Apply clean 22mm² coverslip to slide. Seal edges with rubber cement.
- 7. Co-denature chromosomes/probe on an hotplate at 72°C for 2 minutes, or 83°C for 5 minutes for FFPE.
- 8. Place slide in a sealed humidified slide chamber.
- 9. Incubate at 43°C for 1-2 hours.

Post-Hybridization Washes

- 1. Pre-warm WS1 (0.3% Igepal (Sigma CA-630) or NP-40 / 0.4 x SSC) for 1 hour to 73° C.
- 2. Carefully remove the rubber cement around the coverslip.
- 3. Place slide in WS1 and let stand for exactly 2 minutes.
- 4. Transfer to WS2 at room temperature for 2 minutes.
- 5. Let dry in dark.
- 6. Apply 10 μl DAPI with Antifade and 22mm² coverslip.
- 7. Wait 15-30 minutes then visualize under microscope using the appropriate filter sets.

References

Barch MJ, Knutsen T, Spurbeck JL. The AGT Cytogenetics Laboratory Manual, Third Edition. Lippincott-Raven Philadelphia. 1991.