

Standard Operating Procedure Form		
F-195 RTU FISH Probe Automated Hybridization IFU		
Rev Date: 09/03/2025	Revision: 0	1 of 2
Authorized By: Sofia Badanin	Approved By: Richie Robel	

## SwishProbes<sup>RTU</sup> – Hybridization Protocol

### Notes

- Protocol can be used with all FISH probes – controls, gene specifics, custom FISH probes.
- Solutions can be made prior to the procedure.
- Further optimization of the protocol may be required.

### Required Reagents & Equipment (Not Supplied)

70%, 85%, and 100% EtOH

22x22 mm dimension coverglass

22x50 mm dimension coverglass

Rubber cement

HYBrite or ThermoBrite

Absorbent Material

dH<sub>2</sub>O

Wash Buffer I (SKU: WASH1-004L or WASH1-250M)

Wash Buffer II (SKU: WASH2-004L or WASH2-250M)

DAPI with Antifade (SKU: DAPI-0125)

**Prep:** Pre-soak absorbent material in dH<sub>2</sub>O to place in HYBrite/ThermoBrite. Prepare wash solutions and warm WS1 to 73°C

### Automated Hybridization Procedure:

1. Prepare the program on the HYBrite/ThermoBrite:
  - a. Peripheral Blood/ Bone Marrow Cell Pellets:
    - a. 2 Hour Hybridization:
      - i. Denaturation: 80-83°C for 2-3 minutes; Hybridization: 37°C for 2 hours.
    - b. 16 Hour Hybridization:
      - i. Denaturation: 80-83°C for 2-3 minutes; Hybridization: 37°C for 16 hours.
  - b. FFPE Tissue Samples:
    - a. 2 Hour Hybridization:
      - i. Denaturation: 83°C for 5-7 minutes; Hybridization: 37°C for 2 hours.
    - b. 16 Hour Hybridization:
      - i. Denaturation: 83°C for 5-7 minutes; Hybridization: 37°C for 16 hours.
2. Place pre-soaked absorbent material in HYBrite/ThermoBrite and ensure HYBrite/ThermoBrite is on and warming.
3. Dehydrate slides:
  - a. Place slides in 70% EtOH for 2 minutes.
  - b. Transfer slides in 85% EtOH for 2 minutes.
  - c. Transfer slides in 100% EtOH for 2 minutes.
  - d. Air dry slides.
4. Mix probe thoroughly by vortexing before use.
5. Pipette 10µL of the probe mixture onto the slides over each cellular area and place 22x22 mm coverglass on the slides to cover each cellular area, taking care that there are no air bubbles.
6. Seal the coverglass with rubber cement. Be sure that all edges of coverglass have been sealed with the rubber cement so the probe does not dry out under the coverslip.
7. Place the slides onto the HYBrite/ThermoBrite, close the lid, and begin the program that is specific to your cell type.
8. After incubation period, remove the slides from the HYBrite/ThermoBrite, remove the rubber cement, gently remove the coverglass, and move on to the wash process.

Standard Operating Procedure Form		
F-195 RTU FISH Probe Automated Hybridization IFU		
Rev Date: 09/03/2025	Revision: 0	2 of 2
Authorized By: Sofia Badanin	Approved By: Richie Robel	

9. Wash Slides:

- a. Place slides in Wash Buffer I at 73°C for 2 minutes; agitate slides for the first 10-15 seconds.
- b. Transfer slides into Wash Buffer II at room temperature for 2 minutes; agitate slides for the first 10-15 seconds.
- c. Remove the slides and wipe the backside to dry. Allow the rest of the slide to dry in a dark place.

10. Once the slide is completely dry, apply 10µL of DAPI with Antifade to the cellular area and apply the 22x50 mm coverslip. Store slide in a dark place and wait 15-20 minutes before viewing it under the microscope.