

Standard Operating Procedure Form		
F-83 SwiftFISH Rapid Hybridization Buffer IFU		
Rev Date: 08/11/2025	Revision: 7	1 of 1
Authorized By: Sofia Badanin	Approved By: Richie Robel	

SwiftFISH Rapid Hybridization Buffer

Notes

- The SwiftFISH Rapid Hybridization Buffer decreases hybridization time to 2 hours and is also optimized to perform a traditional 16-hour-hybridization.
- Works for the following sample types: peripheral blood (PB), bone marrow (BM) and formalin-fixed paraffin-embedded (FFPE).
- Works with all FISH probes: controls, gene specific, custom FISH probes.
- Further optimization of the protocol may be required.
- Thaw and mix the buffer well prior to use.

Required Reagents & Equipment (Not Supplied)

Thermobrite

70%, 85%, 100% Ethanol

Wash Solution 1 (WS1) – 0.3% Igepal (Sigma CA-630) or NP-40 / 0.4 x SSC

Wash Solution 2 (WS2) – 0.1% Igepal (Sigma CA-630) or NP-40 / 2 x SSC

DAPI with Antifade

18x18mm, 22x22mm, 22x50mm coverslip

Hybridization Setup

- Take buffer from freezer and bring to room temperature.
- Mix buffer by vortexing for approximately 10 seconds, making sure there is no precipitate.
- Place slides on a warm plate (45°C) for 20 minutes – you can decrease time on the plate by increasing the temperature (i.e. 70°C for 10 minutes).
- Dehydrate the slides in Ethanol series 70%, 85%, and 100%, 2 minutes each. Let dry.
- Add 10 µl probe mixture to slide (2 µl probe + 8 µl SwiftFISH Rapid Hybridization Buffer) using a pipette tip with the end cut if needed as the buffer is thick.
- Apply clean 18mm² or 22mm² coverslip to slide. Seal edges with rubber cement.
- Place slides on Thermobrite/Hybrite with the following setting:
 - Peripheral Blood/ Bone Marrow Cell Pellets:
 - 2 Hour Hybridization:
 - Denaturation: 80-83°C for 2-3 minutes; Hybridization: 37°C for 2 hours.
 - 16 Hour Hybridization:
 - Denaturation: 80-83°C for 2-3 minutes; Hybridization: 37°C for 16 hours.
 - FFPE Tissue Samples:
 - 2 Hour Hybridization:
 - Denaturation: 83°C for 5-7 minutes; Hybridization: 37°C for 2 hours.
 - 16 Hour Hybridization:
 - Denaturation: 83°C for 5-7 minutes; Hybridization: 37°C for 16 hours.

Post-Hybridization Washes

- Pre-warm WS1 (0.3% Igepal (Sigma CA-630) or NP-40 / 0.4 x SSC) for 1 hour to 73°C.
- Carefully remove the rubber cement around the coverslip.
- Place slide in WS1 and let stand for exactly 2 minutes, agitating for the first ~15 seconds.
- Transfer to WS2 at room temperature for 2 minutes, agitating for the first ~15 seconds.
- Let dry in dark.
- Apply 10 µl DAPI with Antifade and 22x50mm coverslip.
- Wait 15-30 minutes then visualize under microscope using the appropriate filter sets.

References

Barch MJ, Knutsen T, Spurbeck JL. The AGT Cytogenetics Laboratory Manual, Third Edition. Lippincott-Raven Philadelphia. 1991.