




The CEP85 FISH probe is designed to hybridize to the CEP85 gene on 1p36.11 and is primarily used for detecting amplifications and deletions associated with the gene. This probe is FISH confirmed on normal peripheral blood metaphase spreads and interphase nuclei. The probe can be labeled in one of five colors, using standard Nick Translation protocols. Each probe is sold in a 20 test kit (approximately 20 slides - 22x22 mm area) and includes hybridization buffer. Please also note that due to design optimizations, prices are subject to change.



| SKU         | Test Kits  | Buffer | Dye Color | Price   |
|-------------|------------|--------|-----------|---------|
| CEP85-20-RE | 20 (40 μL) | 200 μL |           | \$1,500 |
| CEP85-20-OR | 20 (40 μL) | 200 μL |           | \$1,200 |
| CEP85-20-GO | 20 (40 μL) | 200 μL |           | \$1,500 |
| CEP85-20-GR | 20 (40 μL) | 200 μL |           | \$1,500 |
| CEP85-20-AQ | 20 (40 μL) | 200 μL |           | \$1,500 |

## **Ordering Instructions:**

To order the CEP85 FISH Probe, visit https://www.empiregenomics.com/fish-probes/gene/CEP85 or contact our office at (716) 856-3873.