






The FECD2 FISH probe is designed to hybridize to the FECD2 gene on 13pter-q12.13 and is primarily used for detecting amplifications and deletions associated with the gene. This probe is FISH confirmed on normal peripheral blood metaphase spreads and interphase nuclei. The probe can be labeled in one of five colors, using standard Nick Translation protocols. Each probe is sold in a 20 test kit (approximately 20 slides - 22x22 mm area) and includes hybridization buffer. Please also note that due to design optimizations, prices are subject to change.



| SKU | Test Kits | Buffer | Dye Color | Price |
|-------------|------------|--------|---|---------|
| FECD2-20-RE | 20 (40 µL) | 200 µL |  | \$1,500 |
| FECD2-20-OR | 20 (40 µL) | 200 µL |  | \$1,200 |
| FECD2-20-GO | 20 (40 µL) | 200 µL |  | \$1,500 |
| FECD2-20-GR | 20 (40 µL) | 200 µL |  | \$1,500 |
| FECD2-20-AQ | 20 (40 µL) | 200 µL |  | \$1,500 |

Ordering Instructions:

To order the FECD2 FISH Probe, visit <https://www.empiregenomics.com/fish-probes/gene/FECD2> or contact our office at (716) 856-3873.

** For in vitro use only / RUO in US and other countries*