




The NLK FISH probe is designed to hybridize to the NLK gene on 17q11.2 and is primarily used for detecting amplifications and deletions associated with the gene. This probe is FISH confirmed on normal peripheral blood metaphase spreads and interphase nuclei. The probe can be labeled in one of five colors, using standard Nick Translation protocols. Each probe is sold in a 20 test kit (approximately 20 slides - 22x22 mm area) and includes hybridization buffer. Please also note that due to design optimizations, prices are subject to change.



SKU	Test Kits	Buffer	Dye Color	Price
NLK-20-RE	20 (40 μL)	200 μL		\$1,500
NLK-20-OR	20 (40 μL)	200 μL		\$1,200
NLK-20-GO	20 (40 μL)	200 μL		\$1,500
NLK-20-GR	20 (40 μL)	200 μL		\$1,500
NLK-20-AQ	20 (40 μL)	200 μL		\$1,500

## **Ordering Instructions:**

To order the NLK FISH Probe, visit https://www.empiregenomics.com/fish-probes/gene/NLK or contact our office at (716) 856-3873.