



The PDZD2 FISH probe is designed to hybridize to the PDZD2 gene on 5p13.3 and is primarily used for detecting amplifications and deletions associated with the gene. This probe is FISH confirmed on normal peripheral blood metaphase spreads and interphase nuclei. The probe can be labeled in one of five colors, using standard Nick Translation protocols. Each probe is sold in a 20 test kit (approximately 20 slides - 22x22 mm area) and includes hybridization buffer. Please also note that due to design optimizations, prices are subject to change.



SKU	Test Kits	Buffer	Dye Color	Price
PDZD2-20-RE	20 (40 μL)	200 μL		\$1,500
PDZD2-20-OR	20 (40 μL)	200 μL		\$1,200
PDZD2-20-GO	20 (40 μL)	200 μL		\$1,500
PDZD2-20-GR	20 (40 μL)	200 μL		\$1,500
PDZD2-20-AQ	20 (40 μL)	200 μL		\$1,500

## **Ordering Instructions:**

To order the PDZD2 FISH Probe, visit https://www.empiregenomics.com/fish-probes/gene/PDZD2 or contact our office at (716) 856-3873.

\* For in vitro use only | RUO in US and other countries