

The NANOGP7 FISH probe is designed to hybridize to the NANOGP7 gene on 14q32.12 and is primarily used for detecting amplifications and deletions associated with the gene. This probe is FISH confirmed on normal peripheral blood metaphase spreads and interphase nuclei. The probe can be labeled in one of five colors, using standard Nick Translation protocols. Each probe is sold in a 20 test kit (approximately 20 slides - 22x22 mm area) and includes hybridization buffer. Please also note that due to design optimizations, prices are subject to change.

NANOGP7 FISH Probe (14q32.12)









IMAGE NOT DRAWN TO SCALE

Design may change based on optimization, a final design will be provided upon purchase

SKU	Test Kits	Buffer	Dye Color	Price
NANOGP7-20-RE	20 (40 µL)	200 µL		\$1,500
NANOGP7-20-OR	20 (40 µL)	200 µL		\$1,200
NANOGP7-20-GO	20 (40 µL)	200 µL		\$1,500
NANOGP7-20-GR	20 (40 µL)	200 µL		\$1,500
NANOGP7-20-AQ	20 (40 µL)	200 µL		\$1,500

Ordering Instructions:

To order the NANOGP7 FISH Probe, visit <https://www.empiregenomics.com/fish-probes/gene/NANOGP7> or contact our office at (716) 856-3873.